

Aerocrine

NIOX MINO Speaker's video

EPV-000007-00
February 14, 2007

This mpg-file can be used in PowerPoint presentations in order to show a routine examination with NIOX MINO[®] and to demonstrate just how fast and easy it is to use this non-invasive method (F_{ENO} measurement).

There is no voice-over to the mpg-file, but you can hear when the patient is emptying his lungs and the various sounds from the instrument during the procedure. This is of course provided that you are using loudspeakers.

The video is 1:13 min.

Just plug in the disposable filter...touch the display window ...and you are ready to go.

First, ask the patient to empty the lungs and then to make a fast and deep inhalation through the filter to the total lung capacity. There is a NO-scrubber inside the instrument which effectively cleans the air from ambient nitric oxide contamination.

Then, ask the patient to exhale slow and steady through the filter. Light and sound signals as well as an animated symbol on the display help the patient to perform this step correctly.

A "click" sounds... telling you that the measurement was successful. The device will analyze the sample within a few minutes and signal when the result is ready to read.